



- Guaranteed dew point performance
- Design simplicity
- No refrigerants used
- Purge saving & Dew point dependent control system
- Compact digital CO & Dewpoint Monitor/Alarm
- Air quality mandated by NFPA 99
- Bacterial Penetration upto 0.0001%
- LED tower operation display indicating sequence of operation




Trident

**Breathing Air System** 

# **TBAS Medical Dryers**

The TBAS operate with the principle of Pressure Swing Adsorption, hence removing more moisture than refrigeration air dryers, and delivers moisture free dry air consistently irrespective of flow variation.



# **DESCRIPTION OF OPERATION**

### **Drying cycle**

The compressed wet air flows through the pre filter. The Liquid water particles get filtered by the filter. The filtered air flows in to the adsorber tower filled with desiccant where it loses all the moisture to the desiccant. Purified (Moisture and oil free) air further passing through the series filter to achieve the desired quality level as per Breathing Standard.

### **Regeneration Cycle**

In order to remove the moisture during regeneration cycle. Small portion of dry air from the drying tower is passes over the desiccant through the regeneration orifice. This results in complete regeneration of desiccant and ready for the next cycle.

Iten Cod		L mm	B mm	H mm	Flow cfm	End Conn. BSP	Item Code	Model TBAS #Duplex	L mm	B mm	H mm
PM0	10S	625	540	1121	10	1/2"	PM009	10D	1150	800	1168
PM0	20S	625	540	1121	20	1/2"	PM010	20D	1150	800	1168
PM0	<b>3</b> 0S	625	540	1316	30	1/2"	PM011	30D	1150	800	1327
PM0	<b>45</b> S	663	648	1120	45	1/2"	PM012	45D	1150	800	1170
PM0	<b>60</b> S	808	740	1290	60	1"	PM013	60D	1150	838	1327
PM0	100S	823	740	1745	100	1"	PM014	100D	1343	800	1788
PM0	125S	823	740	2055	125	1"	PM015	125D	1343	805	2085
PM0	<b>200</b> S	865	795	1819	200	11/2"	PM016	200D	1369	912	1819
PM04	<b>11</b> 250S	851	797	2074	250	11/2"	PM040	250D	1465	933	2231
PM04	300S	990	992	1947	300	2"	PM039	300D	1777	1060	1947
PM04	<b>375</b> S	990	992	2257	375	2"	PM037	375D	1787	1050	2258

#### Standard Features

- \* Single breathing Air unit with CO and Dew Point Monitor (Simplex)
- # Dual breathing Air unit with CO and Dew Point Monitor (Duplex) One Stand by unit with manual redundancy.

# **Operating Conditions**

Voltage : 100 - 240 VAC, 50/60Hz, 1ph

Power Consumption : 40 Watt Max Maximum working pressure : 16 bar (g) Minimum working pressure : 4 bar (g) Rated inlet Temperature : 38°C Rated inlet pressure : 7 bar(g)

CO Level : 10 PPM (Alarm set point)

as per NFPA 99

Air Quality : 2°C PDP (Alarm set point)

as per NFPA 99

0.01 micron particle level.

Cycle Time : 4 mir Dewpoint dependent change over

The **TBAS** Series has robust valving and includes high quality coalescing pre-filter, a particulate fine-filter, a Activated carbon filter, and Bacterial filter mounted as a standard feature.

### **Coalesing - Pre Filter**

The P-Grade (5 Micron) pre-filter element removes particulates, water aerosols, and oil mist content in the air.

# Fine Filter

The X-Grade (1 Micron) Fine filter efficiently removes fine aerosols and desiccant powder particles from the adsorbent towers.

#### **Activated Carbon Filter**

This filter is used to remove oil and hydrocarbon vapor from the compressed air stream coming from the dryer. Odor and Taste free air at the outlet.

#### **Bacterial Filter**

Bacterial filters (0.01 Micron) provide effective protection against various types of particles including bacteria, viruses, and moisture droplets in the dry air out from the dryer.

These are help to protect the patient, and the breathing circuit, hospital pipe line from contamination, odor and taste free dry air at the outlet.

### **TBAS** controller

TBAS provides complete control of the system with LED indication/description of each step in the sequence of operation. A high quality Carbon Monoxide (CO) Monitor and Dew point monitor ensures the continuous monitoring of CO in PPM level and Water vapor concentrations in dewpoint deg °C in the Digital display to ensure the safety of all users.

## TBAS Outstanding features and applications

## ✓ Alarms

Compact digital carbon monoxide (CO) monitor, Dew point monitor. If the set value exceeds the standard value, it will give the alarm in the front panel.

### ✓ Purge saving

Manual Purge economiser.

Dewpoint sensor for automatic purge control.

### ✓ Construction

Corrosion free aluminium towers and blocks. All pipeline are in SS/Copper.

### ✓ Display

LED indication for tower operation and sequence of operations.

### √ Easy Integration

Easy installation and integration with existing equipment: All system tie-in points are on one side.

# ✓ Safe

Heavy-duty adsorption vessels, designed and certified for an unlimited number of cyclic loads.

### ✓ Customization

An extended list of options allow you to define your specific Trident Twin-Tower PSA breathing air dryer adjusted to your individual need.

The TBAS series breathing air dryer meets/exceeds the breathing air standards below.

Breathing Air Quality:

NFPA99

OSHA : Grade D

European Pharmacopoeia

IS/ISO 7396

Statement of conformity for Dryer and Filter

• ISO 8573 Class - 2 : Outlet air quality

• 97 / 23 / CE : Pressure Directive

• 89 / 392 / CEE : Machinery Directive

73 / 23 / CEE : Low Voltage

CRN

#### **Trident Pneumatics Pvt Ltd**

5/232, K.N.G. Pudur Road, Coimbatore - 641 108, India. Ph: +91-422-2400492, 2401373 Fax: +91-422-2401376

e-mail: sales@tridentpneumatics.com Website: www.tridentpneumatics.com

Andhara Pradesh, Tamilnadu & Maharashtra & After North Region Karnataka Kolkatta Sales Guiarat Telangana M.P Kerala Sales Outlets & Chhattisgarh 4436308 9789480564 9789480564 7397719565 9885445321 9025404538 9867367726 8870010565